

IRAN-LAB.

Iranian Test Machines



HDT/VICAT TESTER



No. 11, Vosoghi St., Metal Industries St., Fath Highway, Tehran

021 66808218

www.iranlab.co

Model	HV-300
Standard number	Vicat : ISO 306 ,ASTM D1525, ISIRI 6982 HDT : ISO 75 ,ASTM D648, ISIRI 6845
Testing method	HDT and VICAT
Temperature range	Room temperature-300°C
Sampeling station	1 up to 3
Deformation Display	Digital with 0.01 mm accuracy
Rate of temperature	2°C/min or 50°C/hr ± 5°C
cooling method	None (recommended by costumer)
Power	220 V 1PH 5 A
Weight	35 Kg