


UNIVERSAL TENSILE MACHINE



Model	TM-500-M
Standard number	ISIRI 6621-1, ISIRI 7276, ISIRI 14471-2 , ASTM D638, ISO 527-1
Maximum force	20 KN (2000 Kg)
Unit (Switchable)	Kg ,N ,N/mm2 ,Kg/cm2
Drive	Electric-motor and inverter
Travel	600 mm by reading accuracy 0.001 mm
Processor	DELTA PLC
Applications	Cellulose, Textile, Polymer, adhesives, Wood
Test Speed	5 - 80 mm/min
Weight (approx.)	80 Kg
Columns	2 rod Ball-screw
Monitoring	DELTA 4.3 inch multi-color touch screen HMI Controller
Column Spacing	350 mm
Grips	1 set tensile wedge grips on machine
Power Supply	220V 5A