

# IRAN-LAB.

Iranian Test Machines



## RUB TESTER



No. 11, Vosoghi St., Metal Industries St., Fath Highway, Tehran

021 66808218

[www.iranlab.co](http://www.iranlab.co)

<b>Model</b>	<b>RT-100</b>
Standard Number	UMA – 5266 , ASTM D5264
Loads	0.9 and 1.8 kgf
Measuring	Magnetic counter
Abrasion length	100 mm
Abraser load	50 * 100 mm
Drive	DC motor