

Melt Flow Indexer (MFI)



Model	MT-300
Standard number	ISO 1133, ASTM D 1238 Method A , ISO 4440 , ISIRI 6980, ISIRI 7175-1
Display	DELTA 4.3 inch Touch screen
Temperature range	Room Temp. ~ 300 `
Load	2160 - 5000~21600 grams
Flow rate	0.15~50 g/10min
Loading Method	Manually
Accessories	Piston rod reamer and cylinder cleaning rod
Auto cut device	Auto/Manual selection
Power	220 V 50 A
Weight	40 KG