

UNIVERSAL TENSILE MACHINE



Model	TM-200-S
Standard number	ISIRI 6621-1, ISIRI 7276 , ASTM D638, ISO 527-1
Monitoring	QCtester 1.0 software on PC
Drive	Servo-motor and servo drive
Applications	Cellulose, Textile, rubber, thin plastic, adhesives and medical
Maximum force	5 KN (500 Kg)
Speed rang	0.05 - 800 mm/min
Processor	Delta PLC Connected to Servo driver
Travel	1200 mm
Device weight	80 Kg
Power	220V 5A