

## Elmendorf Tearing Strength Tester



<b>Model</b>	<b>ET-3200</b>
Standard number	ISIRI 1297, ISO 1974, ASTM D1424 , ISO-6383-2 , ASTM D 1922
Capacity	0~6400 g
Display	DELTA 4.3 inch Touch screen
Test angle	137
Clamping Method	Pneumatic
Distance between clamps	2.8 mm
Grip Distance	Width 25mm, Depth 15.9 mm
PLC	Manufactured by DELTA
Clock pulse	0.001 Second
Result	Angle in degree, Force in gram
Size of specimen	63*80 mm
Length of notch	20 mm
Residual tearing length	43 mm
Weight	30 Kg