

## Temperature and Humidity Chamber



<b>Model</b>	<b>TH-100</b>
Standard Number	ISIRI 2117, ISIRI 106
Temp. range	0°C - 70°C
Temp tolerance	±2.0°C
Temp homogeneity	±3.0°C
Humid Range	30 - 95% R.H
Humid. Tolerance	±5.0% R.H
Humid. Homogeneity	±5.0% R.H
System	Balanced temperature
Heating up time	20°C - 70°C within 30 min
Cooling down time	20°C - 0°C within 40 min
Interior material	Stainless steel or Aluminum
Exterior material	Stainless steel or Painted iron
Insulation	Foam
Refrigeration system	Single stage refrigeration
Power	AC 220 V 8A