

## UVA ACCELERATED WEATHERING TESTER



<b>Model</b>	<b>UFT-340</b>
Standard number	ISO 4892-3
Monitoring	4.3 inch DELTA color touch panel
Programs	Fix and programable
Lamps	8 pcs UVA 340 nm PHILIPS, NARVA
Samples	24 number mounted on aluminum
Alarms	Low water, over heat, finished test time
Heating	PID method by 0.1 degree Accuracy
Body material	Interior tank made by 304 stainless steel, outer side iron cast coated by powder paint