

IRAN-LAB.

Iranian Test Machines



Bursting Strength Tester



No. 11, Vosoghi St., Metal Industries St., Fath Highway, Tehran

021 66808218

www.iranlab.co

Model	BT-2500
Standard number	ISO 2960 , ISIRI 4801, ISIRI 14350
Measuring range	0-5000 Kpa
Accuracy	±0.5% F.S
Speed	170±15 ml/min
Clamp pressure	≥690 Kpa
Diameter of clamp	31.5±0.05
Pneumatic pressure	0.45~0.6 MPa
Dimension	480*420*600 mm
Net Weight	40 kg
Power Supply	220 V 5 A