

## Crush strength tester



<b>Model</b>	<b>CT-200</b>
Standard number	ISIRI 1741, ISIRI 19572-1, ISIRI 18051 , ISIRI 1898 ISIRI 18993, ISO 12192, ISO 3035
Display	DELTA 4.3 inch Touch screen
Driver	Electrical motor powered by gearbox
Capacity	250 KGF
Load accuracy	1 g
Test space	120 mm
Test Speed	12.5 mm/min
Space Between Compression Plate	120 mm
PLC	DELTA 14SS2 PLC
Power Supply	220V 5A